burdinola

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A 3D architectural rendering of a modern building. The structure features a prominent glass and steel framework, with a red brick base. The building has a complex internal layout with multiple levels, including a central atrium and various rooms. The rendering shows the building from an elevated perspective, highlighting its geometric forms and the interplay of light and shadow. The building is situated on a light-colored ground, with some greenery visible in the foreground.

The diagram illustrates a complex air conditioning network. It features a series of air conditioning units, each with specific performance metrics. The units are connected via a network of red and blue pipes, which likely represent different types of refrigerant or air flow. The diagram is overlaid with a grid, and a legend in the top right corner provides additional context.

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Vitrina Gases ELITE 1500 (375 m3/h)

Plano interno

Average concentration in ppm from 60-360s

upper left	upper middle	upper right
0.07	0.05	0.05
0.05	0.01	0.02
lower left	lower middle	lower right
0.01	0.07	< 0.01
<0.01	<0.01	<0.01

Containment factor Cf_i: from 457 to 3200 (limited by detection limit)

Plano externo

Measurement plane: front sash at 500 mm vertical open:

mean concentration in ppm of:			
60s - 360s	361s - 420s	421s - 600s	601s - 780s
<0,01	<0,01	<0,01	0,04
<0,01	<0,01	<0,01	0,02

Containment factor Cf_e: from 1800 to 7200 (limited by detection limit)

Robustez

Measurement plane front sash at 500 mm vertical open:

(Mean concentration in ppm from 60-240s),

	Plane front sash 500 mm vertical open
measurement 1	0,48
measurement 2	0,62
measurement 3	0,38

Containment factor Cf_e: from 116 to 189

NANO Chamber 1800 (470 m3/h)

Plano interno

Average concentration in ppm from 60-360s

upper left	upper middle	upper right
0.027	b.d.l.	b.d.l.
lower left	lower middle	lower right
b.d.l.	b.d.l.	b.d.l.

b.d.l. = below detection limit of the instrument

Plano externo

mean concentration in ppm of:





60s - 360s	361s - 420s	421s - 600s	601s - 660s
0.017	0.024	b.d.l.	b.d.l.

Robustez

Testing of robustness of containment according chapter 5.4

(Mean concentration in ppm from 60-270s),


measurement 1	0.014
measurement 2	0.048
measurement 3	0.015

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BAND	OEB 1	OEB 2	OEB 3	OEB 4	OEB 5	OEB 6	OEB 7
OEL	>1000 - 5000 $\mu\text{g}/\text{m}^3$	>100 - $\leq$ 1000 $\mu\text{g}/\text{m}^3$	>10 - $\leq$ 100 $\mu\text{g}/\text{m}^3$	>1 - $\leq$ 10 $\mu\text{g}/\text{m}^3$	<1.0 - 0.01 $\mu\text{g}/\text{m}^3$	0.01 - 0.001 $\mu\text{g}/\text{m}^3$	<0.001 $\mu\text{g}/\text{m}^3$ - <1 ng/ m^3
Equipment to Use	Ventilation Containment	Ventilation Containment or Flowhoods without downflow (single pass Fume cabinets)	Downflow Booths or VBEs, Flowhoods	VBEs or DFBs with higher containment, Flowhoods with downflow and inflow for small volume	Isolators recommended however if handling less than 3kg and short task duration, low dust cloud potential reverse oRABs possible	Isolators	Isolators
Production OSD				VBEs DFB with high containment screen	WDCI DFB with flexible isolator	WDCI	WDCI
Production common application		Fume Hoods	Downflow Booths		Isolators		
Production Injectable If Aseptic needed (with HPV)				ACTI (Aseptic Containment Isolator) CBI-U (Containment Barrier Isolator Unidirectional) GPPI (General Processing Platform Isolator)	ACTI (Aseptic Containment Isolator) CBI-U (Containment Barrier Isolator Unidirectional) GPPI (General Processing Platform Isolator) with closed transfer	ACTI (Aseptic Containment Isolator) With complete closed transfer (RTPs)	ACTI (Aseptic Containment Isolator) With complete closed transfer (RTPs)
R&D Qc/ IPQC Sampling		FHs (Fume Hoods) EFF (Esco Frontier [®] Floor) EFA (Esco Frontier [®] Acela™) CYT (Cytoculture™ Cytotoxic Safety Cabinets)	FHs (Fume Hoods) EFF (Esco Frontier [®] Floor) EFA (Esco Frontier [®] Acela™) CYT (Cytoculture™ Cytotoxic Safety Cabinets)		GPPI (General Processing Platform Isolator) CBI-T (Containment Barrier Isolator Turbulent) roRABs (Reverse open RABs) VBEs (if small quantity)	CBI-T (Containment Barrier Isolator Turbulent)	CBI-T (Containment Barrier Isolator Turbulent)

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SWAB - SURFACE SAMPLING		RUN 1		VALUE ng / m3
N° TEST		TYOLOGY (ENVIRONMENTAL / PERSONAL)	POSITION	
1C21244 – SW1		Personal	Operator chest	322,4
1C21244 – SW2		Personal	On the right hand	686,4
1C21244 – SW3		Environmental	On the metal surface	20
1C21244 – SW4		Environmental	Middle front of nano chamber	20

TEST 2

AIR SAMPLING POINTS		RUN 2		VALUE ng / m3
N° TEST	S.N. IOM	TYOLOGY (ENVIRONMENTAL / PERSONAL)	POSITION	
1C21244 – SP1	XP2090 - 1	Personal	Operator chest	166,6
1C21244 – SP2	XP2089 - 2	Environmental	On the glass (lateral position)	218,3
1C21244 – SP3	XP2088 - 3	Environmental	Middle front of nano chamber	166,6
1C21244 – SP4	XP2087 - 4	Environmental	Background	166,6

SWAB - SURFACE SAMPLING		RUN 2		VALUE ng / m3
N° TEST		TYOLOGY (ENVIRONMENTAL / PERSONAL)	POSITION	
1C21244 – SW1		Personal	Operator chest	82,1
1C21244 – SW2		Personal	On the right hand	1175,4
1C21244 – SW3		Environmental	On the metal surface	20
1C21244 – SW4		Environmental	Middle front of nano chamber	20

TEST 3

AIR SAMPLING POINTS		RUN 3		VALUE ng / m3
N° TEST	S.N. IOM	TYOLOGY (ENVIRONMENTAL / PERSONAL)	POSITION	
1C21244 – SP1	XP2090 - 1	Personal	Operator chest	1041,6
1C21244 – SP2	XP2089 - 2	Environmental	On the glass (lateral position)	608,3
1C21244 – SP3	XP2088 - 3	Environmental	Middle front of nano chamber	560
1C21244 – SP4	XP2087 - 4	Environmental	Background	381,6

SWAB - SURFACE SAMPLING		RUN 3		VALUE ng / m3
N° TEST		TYOLOGY (ENVIRONMENTAL / PERSONAL)	POSITION	
1C21244 – SW1		Personal	Operator chest	70,6
1C21244 – SW2		Personal	On the right hand	6832
1C21244 – SW3		Environmental	On the metal surface	266
1C21244 – SW4		Environmental	Middle front of nano chamber	51,2

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